

# Bovine BHMT/Betaine--homocysteine S-methyltransferase 1 ELISA Kit

Catalog No: #EK19585

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Bovine BHMT/Betaine--homocysteine S-methyltransferase 1 ELISA Kit          |
| Brief Description  | ELISA Kits   |
| Applications       | ELISA  |
| Species Reactivity | Bovine   |
| Specificity        | Natural and recombinant Bovine Betaine--homocysteine S-methyltransferase 1 |
| Target Name        | Betaine--homocysteine S-methyltransferase 1                                |
| Accession No.      | Q5I597   |
| Uniprot            | Q5I597   |
| GeneID             | 497025;  |

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.86pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BHMT

Note: This product is for in vitro research use only