

# Mouse Cacna1a/Voltage-dependent P/Q-type calcium channel subunit alpha-1A ELISA Kit



Catalog No: #EK19586

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cacna1a/Voltage-dependent P/Q-type calcium channel subunit alpha-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Voltage-dependent P/Q-type calcium channel subunit alpha-1A
Target Name	Voltage-dependent P/Q-type calcium channel subunit alpha-1A
Other Names	Brain calcium channel I, BI, Calcium channel, L type, alpha-1 polypeptide isoform 4, Voltage-gated calcium channel subunit alpha Cav2.1, Caca1a, Cach4, Cacn3, Cacn1a4, Ccha1a
Accession No.	P97445
Uniprot	P97445
GeneID	12286;

## Application Details

Assay Type: Sandwich Detect Range: 7.8-500pg/mL Sensitivity: 2.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cacna1a

Note: This product is for in vitro research use only