

Human ATF1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit



Catalog No: #EK19587

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATF1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cyclic AMP-dependent transcription factor ATF-1
Target Name	Cyclic AMP-dependent transcription factor ATF-1
Other Names	cAMP-dependent transcription factor ATF-1,Activating transcription factor 1,Protein TREB36
Accession No.	P18846
Uniprot	P18846
GeneID	466;

Application Details

Assay Type: Sandwich Detect Range:6.25-400pg/mL Sensitivity:3.125pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ATF1

Note: This product is for in vitro research use only