

Mouse Atf1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit



Catalog No: #EK19588

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Atf1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cyclic AMP-dependent transcription factor ATF-1 |
| Target Name | Cyclic AMP-dependent transcription factor ATF-1 |
| Other Names | cAMP-dependent transcription factor ATF-1,Activating transcription factor 1,TCR-ATF1 |
| Accession No. | P81269 |
| Uniprot | P81269 |
| GeneID | 11908; |

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:5.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Atf1

Note: This product is for in vitro research use only