

Human ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit

Catalog No: #EK19589

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Pituitary adenylate cyclase-activating polypeptide
Target Name	Pituitary adenylate cyclase-activating polypeptide
Other Names	PACAP
Accession No.	P18509
Uniprot	P18509
GeneID	116;

Application Details

Assay Type: Sandwich Detect Range: 6.25-400pg/mL Sensitivity: 3.155pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADCYAP1

Note: This product is for in vitro research use only