

# Mouse Adcyap1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit



Catalog No: #EK19590

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adcyap1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Pituitary adenylate cyclase-activating polypeptide
Target Name	Pituitary adenylate cyclase-activating polypeptide
Other Names	PACAP, Pacap
Accession No.	O70176
Uniprot	O70176
GeneID	11516;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 32pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Adcyap1

Note: This product is for in vitro research use only