

Pig ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit



Catalog No: #EK19591

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pig ADCYAP1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Pituitary adenylate cyclase-activating polypeptide
Target Name	Pituitary adenylate cyclase-activating polypeptide
Other Names	PACAP
Accession No.	P41535
Uniprot	P41535
GeneID	414283;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADCYAP1

Note: This product is for in vitro research use only