

Rat Adcyap1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit



Catalog No: #EK19592

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adcyap1/Pituitary adenylate cyclase-activating polypeptide ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Pituitary adenylate cyclase-activating polypeptide
Target Name	Pituitary adenylate cyclase-activating polypeptide
Other Names	PACAP
Accession No.	P13589
Uniprot	P13589
GeneID	24166;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.11ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adcyap1

Note: This product is for in vitro research use only