

## Human ADCY10/Adenylate cyclase type 10 ELISA Kit

Catalog No: #EK19593

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ADCY10/Adenylate cyclase type 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenylate cyclase type 10
Target Name	Adenylate cyclase type 10
Other Names	AH-related protein, Adenylate cyclase homolog, Germ cell soluble adenylyl cyclase, hsAC, Testicular soluble adenylyl cyclase, SAC
Accession No.	Q96PN6
Uniprot	Q96PN6
GeneID	55811;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.18ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ADCY10

Note: This product is for in vitro research use only