

Human ADCY1/Adenylate cyclase type 1 ELISA Kit

Catalog No: #EK19595

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ADCY1/Adenylate cyclase type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenylate cyclase type 1
Target Name	Adenylate cyclase type 1
Other Names	ATP pyrophosphate-lyase 1, Adenylate cyclase type I, Adenylyl cyclase 1, Ca(2+)/calmodulin-activated adenylyl cyclase
Accession No.	Q08828
Uniprot	Q08828
GeneID	107;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.085ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ADCY1

Note: This product is for in vitro research use only