

## Human AHR/Aryl hydrocarbon receptor ELISA Kit

Catalog No: #EK19597

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human AHR/Aryl hydrocarbon receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aryl hydrocarbon receptor
Target Name	Aryl hydrocarbon receptor
Other Names	Ah receptor, Class E basic helix-loop-helix protein 76, bHLHe76, BHLHE76
Accession No.	P35869
Uniprot	P35869
GeneID	196;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.098ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: AHR

Note: This product is for in vitro research use only