

Rat Ahr/Aryl hydrocarbon receptor ELISA Kit

Catalog No: #EK19598

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ahr/Aryl hydrocarbon receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aryl hydrocarbon receptor
Target Name	Aryl hydrocarbon receptor
Other Names	Ah receptor
Accession No.	P41738
Uniprot	P41738

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.192ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ahr

Note: This product is for in vitro research use only