

Rat Alox5/Arachidonate 5-lipoxygenase ELISA Kit

Catalog No: #EK19599

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Alox5/Arachidonate 5-lipoxygenase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Arachidonate 5-lipoxygenase
Target Name	Arachidonate 5-lipoxygenase
Other Names	5-LO
Accession No.	P12527
Uniprot	P12527
GeneID	25290;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.194ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Alox5

Note: This product is for in vitro research use only