

Mouse Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit



Catalog No: #EK19601

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Arachidonate 15-lipoxygenase B |
| Target Name | Arachidonate 15-lipoxygenase B |
| Other Names | 15-LOX-B, 15-lipoxygenase 2, 15-LOX-2, 8S-lipoxygenase, 8S-LOX, Arachidonate 15-lipoxygenase type II, Alox8 |
| Accession No. | O35936 |
| Uniprot | O35936 |
| GeneID | 11688; |

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Alox15b

Note: This product is for in vitro research use only