

Rat Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit



Catalog No: #EK19602

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Arachidonate 15-lipoxygenase B
Target Name	Arachidonate 15-lipoxygenase B
Other Names	15-LOX-B, 15-lipoxygenase 2, 15-LOX-2, Arachidonate 15-lipoxygenase type II, Linoleate 13-lipoxygenase 15-LOb
Accession No.	Q8K4F2
Uniprot	Q8K4F2
GeneID	266604;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Alox15b

Note: This product is for in vitro research use only