

Human ATM/Serine-protein kinase ATM ELISA Kit

Catalog No: #EK19603

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ATM/Serine-protein kinase ATM ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Serine-protein kinase ATM
Target Name	Serine-protein kinase ATM
Other Names	Ataxia telangiectasia mutated,A-T mutated
Accession No.	Q13315
Uniprot	Q13315
GeneID	472;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.16ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ATM

Note: This product is for in vitro research use only