

Mouse Atm/Serine-protein kinase ATM ELISA Kit

Catalog No: #EK19604

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atm/Serine-protein kinase ATM ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serine-protein kinase ATM
Target Name	Serine-protein kinase ATM
Other Names	Ataxia telangiectasia mutated homolog,A-T mutated homolog
Accession No.	Q62388
Uniprot	Q62388
GeneID	11920;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:40.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Atm

Note: This product is for in vitro research use only