

Human BRCA2/Breast cancer type 2 susceptibility protein ELISA Kit

Catalog No: #EK19605

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BRCA2/Breast cancer type 2 susceptibility protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Breast cancer type 2 susceptibility protein
Target Name	Breast cancer type 2 susceptibility protein
Other Names	Fanconi anemia group D1 protein,FACD,FANCD1
Accession No.	P51587
Uniprot	P51587
GeneID	675;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:40.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:BRCA2

Note: This product is for in vitro research use only