

Mouse Ada/Adenosine deaminase ELISA Kit

Catalog No: #EK19607

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ada/Adenosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenosine deaminase
Target Name	Adenosine deaminase
Other Names	Adenosine aminohydrolase
Accession No.	P03958
Uniprot	P03958
GeneID	11486;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.174ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ada

Note: This product is for in vitro research use only