

## Human CST4/Cystatin-S ELISA Kit

Catalog No: #EK19609

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CST4/Cystatin-S ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cystatin-S
Target Name	Cystatin-S
Other Names	Cystatin-4,Cystatin-SA-III,Salivary acidic protein 1
Accession No.	P01036
Uniprot	P01036
GeneID	1472;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.178ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CST4

Note: This product is for in vitro research use only