

Human CIB1/Calcium and integrin-binding protein 1 ELISA Kit



Catalog No: #EK19610

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CIB1/Calcium and integrin-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calcium and integrin-binding protein 1
Target Name	Calcium and integrin-binding protein 1
Other Names	CIB, Calcium- and integrin-binding protein, CIBP, Calmyrin, DNA-PKcs-interacting protein, Kinase-interacting protein, KIP, SNK-interacting protein 2-28, SIP2-28, CIB, KIP, PRKDCIP
Accession No.	Q99828
Uniprot	Q99828
GeneID	10519;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.47ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CIB1

Note: This product is for in vitro research use only