

# Mouse Cib1/Calcium and integrin-binding protein 1 ELISA Kit

Catalog No: #EK19611

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mouse Cib1/Calcium and integrin-binding protein 1 ELISA Kit                                     |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Ms  |
| Specificity        | Natural and recombinant Mouse Calcium and integrin-binding protein 1                            |
| Target Name        | Calcium and integrin-binding protein 1  |
| Other Names        | CIB, Calmyrin, DNA-PKcs-interacting protein, Kinase-interacting protein, KIP, Cib, Kip, Prkdcip |
| Accession No.      | Q9Z0F4  |
| Uniprot            | Q9Z0F4  |
| GeneID             | 23991;  |

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.435ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cib1

Note: This product is for in vitro research use only