

Rat Cib1/Calcium and integrin-binding protein 1 ELISA Kit



Catalog No: #EK19612

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cib1/Calcium and integrin-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calcium and integrin-binding protein 1
Target Name	Calcium and integrin-binding protein 1
Other Names	CIB, Calmyrin, DNA-PKcs-interacting protein, Kinase-interacting protein, KIP, Cib, Kip, Prkdcip
Accession No.	Q9R010
Uniprot	Q9R010
GeneID	81823;

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.267 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cib1

Note: This product is for in vitro research use only