

## Human CSN3/Kappa-casein ELISA Kit

Catalog No: #EK19613

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CSN3/Kappa-casein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Kappa-casein
Target Name	Kappa-casein
Other Names	CASK,CSN10,CSNK
Accession No.	P07498
Uniprot	P07498
GeneID	1448;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.066ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CSN3

Note: This product is for in vitro research use only