

## Mouse Cmpk1/UMP-CMP kinase ELISA Kit

Catalog No: #EK19619

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cmpk1/UMP-CMP kinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse UMP-CMP kinase
Target Name	UMP-CMP kinase
Other Names	Deoxycytidylate kinase,CK,Nucleoside-diphosphate kinase,Uridine monophosphate/cytidine monophosphate kinase,UMP/CMP kinase,Cmk,Cmpk,Uck,Umk,Umpk
Accession No.	Q9DBP5
Uniprot	Q9DBP5
GeneID	66588;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cmpk1

Note: This product is for in vitro research use only