

Human CGREF1/Cell growth regulator with EF hand domain protein 1 ELISA Kit



Catalog No: #EK19620

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CGREF1/Cell growth regulator with EF hand domain protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cell growth regulator with EF hand domain protein 1
Target Name	Cell growth regulator with EF hand domain protein 1
Other Names	Cell growth regulatory gene 11 protein,Hydrophobestin,CGR11
Accession No.	Q99674
Uniprot	Q99674
GeneID	10669;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.194ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CGREF1

Note: This product is for in vitro research use only