

Human BCAP31/B-cell receptor-associated protein 31 ELISA Kit



Catalog No: #EK19621

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BCAP31/B-cell receptor-associated protein 31 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human B-cell receptor-associated protein 31
Target Name	B-cell receptor-associated protein 31
Other Names	BCR-associated protein 31,6C6-AG tumor-associated antigen,Protein CDM,p28,BAP31,DXS1357E
Accession No.	P51572
Uniprot	P51572
GeneID	10134;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:32pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BCAP31

Note: This product is for in vitro research use only