

Mouse Bcap31/B-cell receptor-associated protein 31 ELISA Kit



Catalog No: #EK19622

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bcap31/B-cell receptor-associated protein 31 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell receptor-associated protein 31
Target Name	B-cell receptor-associated protein 31
Other Names	BCR-associated protein 31,p28,Bap31
Accession No.	Q61335
Uniprot	Q61335
GeneID	27061;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bcap31

Note: This product is for in vitro research use only