

## MYO1F Antibody

Catalog No: #44670

Package Size: #44670-1 50ul #44670-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MYO1F Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYO1F antibody detects endogenous levels of total MYO1F
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYO1F
Target Name	MYO1F
Other Names	Myo1f antibody MYO1F_HUMAN antibody Myosin ID antibody Myosin IF antibody Myosin-Ie antibody Unconventional myosin-Ie antibody
Accession No.	Swiss-Prot#: O00160 NCBI Gene ID:4542
Uniprot	O00160
GeneID	4542;
Calculated MW	125 kDa.
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only