

Human CST9/Cystatin-9 ELISA Kit

Catalog No: #EK19624

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Human CST9/Cystatin-9 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Cystatin-9 |
| Target Name | Cystatin-9 |
| Other Names | Cystatin-like molecule,CLM,CTES7A |
| Accession No. | Q5W186 |
| Uniprot | Q5W186 |
| GeneID | 128822; |

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CST9

Note: This product is for in vitro research use only