

## Mouse Cst9/Cystatin-9 ELISA Kit

Catalog No: #EK19625

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cst9/Cystatin-9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cystatin-9
Target Name	Cystatin-9
Other Names	Testatin, Cresp
Accession No.	Q9Z0H6
Uniprot	Q9Z0H6
GeneID	13013;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.162ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cst9

Note: This product is for in vitro research use only