

Mouse Asrgl1/Isoaspartyl peptidase/L-asparaginase ELISA Kit



Catalog No: #EK19627

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Asrgl1/Isoaspartyl peptidase/L-asparaginase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Isoaspartyl peptidase/L-asparaginase
Target Name	Isoaspartyl peptidase/L-asparaginase
Other Names	Asparaginase-like protein 1, Beta-aspartyl-peptidase, Isoaspartyl dipeptidase, L-asparagine amidohydrolase
Accession No.	Q8C0M9
Uniprot	Q8C0M9
GeneID	66514;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.057ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Asrgl1

Note: This product is for in vitro research use only