

# Mouse Apeh/Acylamino-acid-releasing enzyme ELISA Kit



Catalog No: #EK19628

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Apeh/Acylamino-acid-releasing enzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Acylamino-acid-releasing enzyme
Target Name	Acylamino-acid-releasing enzyme
Other Names	AARE,Acyl-peptide hydrolase,APH,Acylaminoacyl-peptidase
Accession No.	Q8R146
Uniprot	Q8R146
GeneID	235606;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.097ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Apeh

Note: This product is for in vitro research use only