

Mouse Ambn/Ameloblastin ELISA Kit

Catalog No: #EK19629

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ambn/Ameloblastin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Ameloblastin
Target Name	Ameloblastin
Accession No.	O55189
Uniprot	O55189
GeneID	11698;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.068ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ambn

Note: This product is for in vitro research use only