

# Mouse Akt1s1/Proline-rich AKT1 substrate 1 ELISA Kit

Catalog No: #EK19631

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Akt1s1/Proline-rich AKT1 substrate 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Proline-rich AKT1 substrate 1
Target Name	Proline-rich AKT1 substrate 1
Other Names	Proline-rich AKT substrate,Pras
Accession No.	Q9D1F4
Uniprot	Q9D1F4
GeneID	67605;

## Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.66pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Akt1s1

Note: This product is for in vitro research use only