

Rat Aktip/AKT-interacting protein ELISA Kit

Catalog No: #EK19632

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Aktip/AKT-interacting protein ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat AKT-interacting protein |
| Target Name | AKT-interacting protein |
| Other Names | Fused toes protein homolog,Fts |
| Accession No. | Q5FVH4 |
| Uniprot | Q5FVH4 |
| GeneID | 291906; |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aktip

Note: This product is for in vitro research use only