

## UNC5C Antibody

Catalog No: #44671

Package Size: #44671-1 50ul #44671-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	UNC5C Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	UNC5C antibody detects endogenous levels of total UNC5C
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human UNC5C
Target Name	UNC5C
Other Names	homolog of C. elegans transmembrane receptor Unc5 antibody Netrin receptor UNC5C antibody Protein unc 5 homolog C antibody Protein unc-5 homolog 3 antibody Protein unc-5 homolog C antibody Unc 5 homolog 3 antibody Unc 5 homolog C antibody Unc5 (C.elegans homolog) c antibody Unc5c antibody UNC5C_HUMA
Accession No.	Swiss-Prot#: O95185NCBI Gene ID:8633
Uniprot	O95185
GeneID	8633;
Calculated MW	103 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only