

Human AICDA/Single-stranded DNA cytosine deaminase ELISA Kit

Catalog No: #EK19633

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AICDA/Single-stranded DNA cytosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Single-stranded DNA cytosine deaminase
Target Name	Single-stranded DNA cytosine deaminase
Other Names	Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase,AID
Accession No.	Q9GZX7
Uniprot	Q9GZX7
GeneID	57379;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.05ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AICDA

Note: This product is for in vitro research use only