

Mouse Aicda/Single-stranded DNA cytosine deaminase ELISA Kit

Catalog No: #EK19634

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aicda/Single-stranded DNA cytosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Single-stranded DNA cytosine deaminase
Target Name	Single-stranded DNA cytosine deaminase
Other Names	Activation-induced cytidine deaminase,AID,Cytidine aminohydrolase,Aid
Accession No.	Q9WVE0
Uniprot	Q9WVE0
GeneID	11628;

Application Details

Assay Type: Sandwich Detect Range:1.56-100U/L Sensitivity:0.73U/L Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aicda

Note: This product is for in vitro research use only