

Human AGPAT1/1-acyl-sn-glycerol-3-phosphate acyltransferase alpha ELISA Kit



Catalog No: #EK19635

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human AGPAT1/1-acyl-sn-glycerol-3-phosphate acyltransferase alpha ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human 1-acyl-sn-glycerol-3-phosphate acyltransferase alpha
Target Name	1-acyl-sn-glycerol-3-phosphate acyltransferase alpha
Other Names	1-acylglycerol-3-phosphate O-acyltransferase 1,1-AGP acyltransferase 1,Lysophosphatidic acid acyltransferase alpha,LPAAT-alpha,Protein G15,G15
Accession No.	Q99943
Uniprot	Q99943
GeneID	10554;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000U/L Sensitivity:7.2U/L Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AGPAT1

Note: This product is for in vitro research use only