

## Mouse Agmat/Agmatinase, mitochondrial ELISA Kit

Catalog No: #EK19636

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Agmat/Agmatinase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Agmatinase, mitochondrial
Target Name	Agmatinase, mitochondrial
Other Names	Agmatine ureohydrolase,AUH
Accession No.	A2AS89
Uniprot	A2AS89
GeneID	75986;

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Agmat

Note: This product is for in vitro research use only