

# Mouse Adamts8/A disintegrin and metalloproteinase with thrombospondin motifs 8 ELISA Kit

Catalog No: #EK19639

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adamts8/A disintegrin and metalloproteinase with thrombospondin motifs 8 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse A disintegrin and metalloproteinase with thrombospondin motifs 8
Target Name	A disintegrin and metalloproteinase with thrombospondin motifs 8
Other Names	ADAM-TS 8,METH-2
Accession No.	P57110
Uniprot	P57110

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.054ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Adamts8

Note: This product is for in vitro research use only