

## Human CA14/Carbonic anhydrase 14 ELISA Kit

Catalog No: #EK19640

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CA14/Carbonic anhydrase 14 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carbonic anhydrase 14
Target Name	Carbonic anhydrase 14
Other Names	Carbonate dehydratase XIV, Carbonic anhydrase XIV, CA-XIV, UNQ690/PRO1335
Accession No.	Q9ULX7
Uniprot	Q9ULX7
GeneID	23632;

## Application Details

Assay Type: Sandwich Detect Range: 12.5-800pg/mL Sensitivity: 6.25pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CA14

Note: This product is for in vitro research use only