

Mouse Ca14/Carbonic anhydrase 14 ELISA Kit

Catalog No: #EK19641

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ca14/Carbonic anhydrase 14 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 14
Target Name	Carbonic anhydrase 14
Other Names	Carbonate dehydratase XIV, Carbonic anhydrase XIV, CA-XIV, Car14, Catm
Accession No.	Q9WVT6
Uniprot	Q9WVT6
GeneID	23831;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.182ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ca14

Note: This product is for in vitro research use only