

Mouse Cited2/Cbp/p300-interacting transactivator 2 ELISA Kit



Catalog No: #EK19642

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cited2/Cbp/p300-interacting transactivator 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cbp/p300-interacting transactivator 2
Target Name	Cbp/p300-interacting transactivator 2
Other Names	MSG-related protein 1,MRG-1,P35srj,Mrg1,Msg2
Accession No.	O35740
Uniprot	O35740
GeneID	17684;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.052ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cited2

Note: This product is for in vitro research use only