

Mouse Ca9/Carbonic anhydrase 9 ELISA Kit

Catalog No: #EK19643

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ca9/Carbonic anhydrase 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 9
Target Name	Carbonic anhydrase 9
Other Names	Carbonate dehydratase IX,Carbonic anhydrase IX,CA-IX,Membrane antigen MN homolog,Car9
Accession No.	Q8VHB5
Uniprot	Q8VHB5
GeneID	230099;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.14ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ca9

Note: This product is for in vitro research use only