

## Rat Chat/Choline O-acetyltransferase ELISA Kit

Catalog No: #EK19644

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Chat/Choline O-acetyltransferase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Choline O-acetyltransferase
Target Name	Choline O-acetyltransferase
Other Names	CHOACTase
Accession No.	P32738
Uniprot	P32738

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.438ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Chat

Note: This product is for in vitro research use only