

Rat Arhgef7/Rho guanine nucleotide exchange factor 7 ELISA Kit



Catalog No: #EK19645

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Arhgef7/Rho guanine nucleotide exchange factor 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Rho guanine nucleotide exchange factor 7
Target Name	Rho guanine nucleotide exchange factor 7
Other Names	Beta-Pix,PAK-interacting exchange factor beta,Pak3bp, Pixb
Accession No.	O55043
Uniprot	O55043
GeneID	114559;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.23ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Arhgef7

Note: This product is for in vitro research use only