

Human BAG6/Large proline-rich protein BAG6 ELISA Kit



Catalog No: #EK19646

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human BAG6/Large proline-rich protein BAG6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Large proline-rich protein BAG6
Target Name	Large proline-rich protein BAG6
Other Names	BAG family molecular chaperone regulator 6,BCL2-associated athanogene 6,BAG-6,HLA-B-associated transcript 3,Protein G3,Protein Scythe,BAT3,G3
Accession No.	P46379
Uniprot	P46379
GeneID	7917;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.43ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:BAG6

Note: This product is for in vitro research use only