

## Mouse Carns1/Carnosine synthase 1 ELISA Kit

Catalog No: #EK19647

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Carns1/Carnosine synthase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carnosine synthase 1
Target Name	Carnosine synthase 1
Other Names	ATP-grasp domain-containing protein 1, Atpgd1, Kiaa1394
Accession No.	Q6ZPS2
Uniprot	Q6ZPS2
GeneID	107239;

## Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.67ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Carns1

Note: This product is for in vitro research use only